



- Automated lab with many available kits



- used with consumable immuSpeed Chips



- driven by easy-to-use immuSpeed Soft



immuSpeed™

Automated:

- mixing
- dilution
- incubation
- washing
- enzymatic readout
- channel regeneration

Your automated ELISA partner

microfluidics for rapid analysis

Features:

- **Fully automated**
Quantitative 96 test results in 1 hour, including incubation
- **Efficient**
Reagent and bead consumption is reduced to a minimum, sample requirement is a mere 10 µl
- **Open**
System runs with conventional reagents and commercially available kits
- **Flexible**
Various assays can be run simultaneously and in different (i.e. urine, plasma, CSF) matrixes

Diagnostic kits

in preparation:

TSH
T3
T4
Testosterone
Progesterone
Prolactin
LH
FSH
hCG
AFP

diagnoSwiss™

*For additional information, please contact:
info@diagnoswiss.com or +41 24 471 4900*